

Venous Return Cannulae, Dual Stage, Low Profile, Wire-reinforced  
Tubing, 32/40 Fr, LRD-61140



The Dual Stage Venous Return Cannula is indicated for single tube venous drainage from the right atrium and vena cava during cardiopulmonary bypass surgery. The Dual Stage Cannula has a proximal flexible polyvinyl chloride tubing segment for drainage of the right atrium. The catheter body tapers to a flexible polyvinyl chloride tubing distal segment that ends in a straight fluted distal tip for drainage of the vena cava. The proximal and distal segments of the catheter are available in wire or non-wire reinforced tubing.

ORDER CODE	LRD-61140
SIZE	32/40 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	41
QTY PER BOX	10
LATEX-FREE	yes