

Venous Return Cannulae, Dual Stage, Wire-reinforced Tubing, 36/50 Fr, RDS-61150



The Dual Stage Venous Return Cannula is indicated for single tube venous drainage from the right atrium and vena cava during cardiopulmonary bypass surgery. The Dual Stage Cannula has a proximal flexible polyvinyl chloride tubing segment for drainage of the right atrium. The catheter body tapers to a flexible polyvinyl chloride tubing distal segment that ends in a straight fluted distal tip for drainage of the vena cava. The proximal and distal segments of the catheter are available in wire or non-wire reinforced tubing.

ORDER CODE	RDS-61150
SIZE	36/50 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	41
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes