

# Venous Return Cannulae, Single Stage, Right Angle Lighthouse Tip, Wire-reinforced Tubing, 32 Fr, RV-41032



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

Order Code	RV-41032
------------	----------

Size	32 Fr
------	-------

Available Coated	no
------------------	----

Qty Per Box	10
-------------	----

Latex-Free	yes
------------	-----

Connector Acceptance	3/8 in
----------------------	--------