

Venous Return Cannulae, Triple Stage, 29/29/29 Fr, RTS-11129



The One Piece Tri Stage Venous Return Cannula is indicated for single tube venous drainage from the right atrium and vena cava during cardiopulmonary bypass surgery. The cannulae body tapers to a flexible polyvinyl chloride tubing distal segment that ends in a straight fluted distal tip for drainage of the vena cava

ORDER CODE	RTS-11129
SIZE	29/29/29 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	41
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes