

# Venous Return Cannulae, Triple Stage, 29/37/37 Fr, RTS-13129



The One Piece Tri Stage Venous Return Cannula is indicated for single tube venous drainage from the right atrium and vena cava during cardiopulmonary bypass surgery. The cannulae body tapers to a flexible polyvinyl chloride tubing distal segment that ends in a straight fluted distal tip for drainage of the vena cava

<b>ORDER CODE</b>	RTS-13129
<b>SIZE</b>	29/37/37 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	41
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes