

# Advanced Circulatory Support, Mechanical Circulatory Support, ProtekDuo, 16/29Fr, 5140-4629



The ProtekDuo dual lumen cannula provides the flexibility to do what's right for each unique patient. It is inserted percutaneously at the right internal jugular vein and it is uniquely positioned to drain deoxygenated blood from the right atrium and return it via the pulmonary artery.

<b>ORDER CODE</b>	5140-4629
<b>SIZE</b>	16/29Fr
<b>AVAILABLE COATED</b>	No
<b>EFFECTIVE LENGTH (CM)</b>	28 proximal / 46 distal
<b>QTY PER BOX</b>	1
<b>LATEX-FREE</b>	yes