

Advanced Circulatory Support, Mechanical Circulatory Support, ProtekSolo Transseptal, 21 Fr, 5140-6217



The ProtekSolo Transseptal cannula provides a powerful cannulation strategy for unloading the left ventricle. The ProtekSolo Transseptal is inserted at the femoral vein and placed in the left atrium via a transseptal procedure; oxygenated blood is returned via the femoral artery.

ORDER CODE	5140-6217
SIZE	21 Fr
AVAILABLE COATED	No
EFFECTIVE LENGTH (CM)	72
QTY PER BOX	1
LATEX-FREE	yes