

# Advanced Circulatory Support, Mechanical Circulatory Support, ProtekSolo Transseptal, 21 Fr, 5140-6221



The ProtekSolo Transseptal cannula provides a powerful cannulation strategy for unloading the left ventricle. The ProtekSolo Transseptal is inserted at the femoral vein and placed in the left atrium via a transseptal procedure; oxygenated blood is returned via the femoral artery.

<b>ORDER CODE</b>	5140-6221
<b>SIZE</b>	21 Fr
<b>AVAILABLE COATED</b>	No
<b>EFFECTIVE LENGTH (CM)</b>	62
<b>QTY PER BOX</b>	1
<b>LATEX-FREE</b>	yes