

Arterial Cannulae, 135 Curved Tip, with Suture Collar and Suture Ring ,
Wire-reinforced Tubing, 12 Fr/3.8 mm, A900-155



The Pediatric cannulae are intended to be used in both the venous and the arterial line of the extracorporeal circuit during cardiopulmonary surgery for periods of up to six hours. Specifically, Arterial Cannulae are designed to return blood to the patient during cardiopulmonary bypass procedures. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

ORDER CODE	A900-155
SIZE	12 Fr/3.8 mm
AVAILABLE COATED	Depending on country registration
EFFECTIVE LENGTH (CM)	22
QTY PER BOX	10
LATEX-FREE	yes