

# Suction Products, Suction Wands, Rigid Fluted Tip with Metal Tube and Plastic Handle, 12.8 Fr, SU- 10101



Suction Wands and Suction Sumps are intended to be used to remove fluids from the surgical field. All of these devices are intended to be used for 6 hours or less. The devices are indicated for use in adult and pediatric patients undergoing surgical procedures requiring cardiopulmonary bypass.

Order Code SU-10101

Size 12.8 Fr

Available Coated no

Qty Per Box 10

Latex-Free yes