

Venous Return Cannulae, Single Stage, Right Angle Lighthouse Tip,
Wire-reinforced Tubing, 24 Fr, RV-41024



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

| | |
|----------------------|----------|
| ORDER CODE | RV-41024 |
| SIZE | 24 Fr |
| AVAILABLE COATED | no |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |
| CONNECTOR ACCEPTANCE | 1/4"-3/8 |