

# Venous Return Cannulae, Single Stage, Straight Lighthouse Tip, Wire-reinforced Tubing, 24 Fr, RV-40024



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

<b>ORDER CODE</b>	RV-40024
<b>SIZE</b>	24 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	36
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>CONNECTOR ACCEPTANCE</b>	1/4" - 3/8