

# Vents, PVC Straight, Bullet Tip, with Guidewire Stylet, 10 Fr, VT-84410



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

<b>ORDER CODE</b>	VT-84410
<b>SIZE</b>	10 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	31,8
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes