


Vents, PVC Right Angle, Bullet Tip, without Stylet, 10 Fr, VT-84510



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

|                       |          |
|-----------------------|----------|
| ORDER CODE            | VT-84510 |
| SIZE                  | 10 Fr    |
| AVAILABLE COATED      | no       |
| EFFECTIVE LENGTH (CM) | 31,8     |
| QTY PER BOX           | 10       |
| LATEX-FREE            | yes      |
| NO DOP                | yes      |