

Arterial Cannulae, Flexible, Curved Tip with Suture Ring, Non-reinforced
Tubing, 18 Fr/ 6.0 mm, NA-5526



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

ORDER CODE	NA-5526
SIZE	18 Fr/ 6.0 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	25
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us