

Arterial Cannulae, Flexible, Straight Tip with Moveable Suture Ring, Non-reinforced Tubing, 24 Fr/ 8.0 mm, NA-3418



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

ORDER CODE	NA-3418
SIZE	24 Fr/ 8.0 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	25
QTY PER BOX	10
LATEX-FREE	yes



**Contact us** Complete the form to get in touch with a LivaNova representative.

Contact Us