

Arterial Cannulae, Flexible, Straight Tip with Moveable Suture Ring, Non-reinforced Tubing, 18 Fr/ 6.0 mm, NA-3436



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

| | |
|------------------------------|---------------|
| ORDER CODE | NA-3436 |
| SIZE | 18 Fr/ 6.0 mm |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 25 |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |