

Arterial Cannulae, Flexible, Straight Tip with Moveable Suture Ring,  
Non-reinforced Tubing, 18 Fr/ 6.0 mm, NA-3436



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

<b>ORDER CODE</b>	NA-3436
<b>SIZE</b>	18 Fr/ 6.0 mm
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	25
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes