

Arterial Cannulae, Flexible, Straight Tip with Suture Ring, Non-reinforced
Tubing, 21 Fr/ 7.0 mm, NA-4527



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

ORDER CODE	NA-4527
SIZE	21 Fr/ 7.0 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	25
QTY PER BOX	10
LATEX-FREE	yes