

# Arterial Cannulae, Optiflow Arterial, Curved Tip, Wire- reinforced Tubing, 21 Fr/ 7.0 mm, A282-70B



The arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary bypass surgery for periods of up to six hours.

Order Code                    A282-70B

Size                             21 Fr/ 7.0 mm

Available Coated            Depending on country registration

Effective Length (cm) 24

Qty Per Box 10

Latex-Free yes

No DOP yes