

# Arterial Cannulae, Optiflow Arterial, Straight Tip, Wire- reinforced Tubing, 24 Fr/ 8.0 mm, A292-80C



The arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary bypass surgery for periods of up to six hours.

Order Code                      A292-80C

Size                                24 Fr/ 8.0 mm

Available Coated                Depending on country registration

Effective Length (cm) 25

Qty Per Box 10

Latex-Free yes

No DOP yes