

Arterial Cannulae, PureFlex, Curved Tip, With Suture Collar, Wire-reinforced Tubing, 18 Fr/6.0 mm, RA18-1BC



The arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary bypass surgery for periods of up to six hours.

Order Code	RA18-1BC
Size	18 Fr/6.0 mm
Available Coated	no
Effective Length (cm)	21.6
Qty Per Box	10
Latex-Free	yes