

# Arterial Cannulae, PureFlex, Curved Tip, With Suture Collar, Wire-reinforced Tubing, 20 Fr/6.5 mm, RA20-1BB



The arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary bypass surgery for periods of up to six hours.

|                       |              |
|-----------------------|--------------|
| Order Code            | RA20-1BB     |
| Size                  | 20 Fr/6.7 mm |
| Available Coated      | no           |
| Effective Length (cm) | 21.6         |
| Qty Per Box           | 10           |
| Latex-Free            | yes          |