

Arterial Cannulae, PureFlex, Curved Tip, With Suture Collar, Wire-reinforced Tubing, 22 Fr/ 7.3 mm, RA22-1BC



The Arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary surgery for periods of up to six hours. The Pureflex Arterial cannulae consist of flexible tubing (tapered to a curved or straight flow tip) made of PVC with wire reinforcement

ORDER CODE	RA22-1BC
SIZE	22 Fr/ 7.3 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	21.6
QTY PER BOX	10
LATEX-FREE	yes



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us