

Arterial Cannulae, PureFlex, Curved Tip, With Suture Flange, Wire-reinforced Tubing, 18 Fr/6.0 mm, RA18-2BB



The Arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary surgery for periods of up to six hours. The Pureflex Arterial cannulae consist of flexible tubing (tapered to a curved or straight flow tip) made of PVC with wire reinforcement

ORDER CODE	RA18-2BB
SIZE	18 Fr/6.0 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	21.6
QTY PER BOX	10
LATEX-FREE	yes



**Contact us** Complete the form to get in touch with a LivaNova representative.

Contact Us