

Cardioplegia Cannulae, Antegrade Aortic Root, with and without Vent Line, 12 Ga (9 Fr), AR-11112



The Aortic Root Cannula is designed to deliver cardioplegia solution during cardiopulmonary bypass surgery and aspirate air from the aorta at the end of the cardiopulmonary bypass procedure. The Aortic Root Cannula with vent line is designed to deliver cardioplegia solution and vent the heart during cardiopulmonary bypass surgery.

ORDER CODE	AR-11112
SIZE	12 Ga (9 Fr)
INSERTION NEEDLE	14 Ga - Green
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	15
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes