

## Cardioplegia Cannulae, Antegrade Aortic Root, with and without Vent Line, 16 Ga (5 Fr), AR-11116



The Aortic Root Cannula is designed to deliver cardioplegia solution during cardiopulmonary bypass surgery and aspirate air from the aorta at the end of the cardiopulmonary bypass procedure. The Aortic Root Cannula with vent line is designed to deliver cardioplegia solution and vent the heart during cardiopulmonary bypass surgery.

|                              |              |
|------------------------------|--------------|
| <b>ORDER CODE</b>            | AR-11116     |
| <b>SIZE</b>                  | 16 Ga (5 Fr) |
| <b>INSERTION NEEDLE</b>      | 18 Ga - Pink |
| <b>AVAILABLE COATED</b>      | no           |
| <b>EFFECTIVE LENGTH (CM)</b> | 15           |
| <b>QTY PER BOX</b>           | 10           |
| <b>LATEX-FREE</b>            | yes          |
| <b>NO DOP</b>                | yes          |