

Cardioplegia Cannulae, Antegrade Aortic Root, with Long Vent Line, 14 Ga (7 Fr), AR-13114



Livanova Antegrade Cardioplegia Aortic Root Cannulae are intended for rapid and secure cardioplegia perfusion of the coronary arteries. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

ORDER CODE	AR-13114
SIZE	14 Ga (7 Fr)
INSERTION NEEDLE	16 Ga - White
EFFECTIVE LENGTH (CM)	15
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes