

# Cardioplegia Cannulae, Antegrade Aortic Root, with Vent Line and Two Colored Clamps, 14 Ga (7 Fr), AR-19114



Livanova Antegrade Cardioplegia Aortic Root Cannulae are intended for rapid and secure cardioplegia perfusion of the coronary arteries. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

<b>ORDER CODE</b>	AR-19114
<b>SIZE</b>	14 Ga (7 Fr)
<b>INSERTION NEEDLE</b>	16 Ga - White
<b>EFFECTIVE LENGTH (CM)</b>	15
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes