

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Conical Tip, Stainless Steel Shaft, 1/4"" Tubing Connector, 11 Fr, P608-35



The Coronary perfusion cannulae are intended to individually cannulate the coronary ostia and to selectively deliver cardioplegic solution during cardiopulmonary bypass surgery for periods up to six hours. Depending on country registration, an optional phospholipid surface treatment may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative

|                  |                                   |
|------------------|-----------------------------------|
| ORDER CODE       | P608-35                           |
| SIZE             | 11 Fr                             |
| AVAILABLE COATED | Depending on country registration |
| QTY PER BOX      | 10                                |
| LATEX-FREE       | yes                               |
| NO DOP           | yes                               |