

## Cardioplegia Cannulae, Flexible Vessels, VC-12000



Livanova Vessel Cannulae are intended for antegrade selective vein graft. They are made of a soft PVC tube, with a soft and flexible distal tip, designed to perfuse vessels from 3.0 to 3.6 mm in diameter. They are available in blunt tip configuration, with or without an integral one-way-check silicone valve.

<b>ORDER CODE</b>	VC-12000
<b>EFFECTIVE LENGTH (CM)</b>	5,1
<b>QTY PER BOX</b>	20
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes