

# Cardioplegia Cannulae, Flexible Vessels, VC-12000



Livanova Vessel Cannulae are intended for antegrade selective vein graft. They are made of a soft PVC tube, with a soft and flexible distal tip, designed to perfuse vessels from 3.0 to 3.6 mm in diameter. They are available in blunt tip configuration, with or without an integral one-way-check silicone valve.

Order Code	VC-12000
------------	----------

Effective Length (cm)	5.1
-----------------------	-----

Qty Per Box	20
-------------	----

Latex-Free	yes
------------	-----

No DOP	yes
--------	-----