

Cardioplegia Cannulae, Retrograde, Manually-Inflating Silicone Balloon, 15 Fr, RCM-14115



The Retrograde Cardioplegia Cannulae are indicated for the infusion of blood or cardioplegia solution into the coronary venous system, intraoperatively. They are available with either a self- or manually-inflating balloon

ORDER CODE	RCM-14115
SIZE	15 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	30,5
QTY PER BOX	10
LATEX-FREE	yes
CONNECTOR	Luer lock