

# Cardioplegia Cannulae, Retrograde, Manually-Inflating Silicone Balloon, 15 Fr, RCM-14215



The Retrograde Cardioplegia Cannulae are indicated for the infusion of blood or cardioplegia solution into the coronary venous system, intraoperatively. They are available with either a self- or manually-inflating balloon

<b>ORDER CODE</b>	RCM-14215
<b>SIZE</b>	15 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	30,5
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>CONNECTOR</b>	Luer lock