

Cardioplegia Cannulae, Retrograde, Manually-Inflating Silicone Balloon, 15 Fr, RCM-14315



The Retrograde Cardioplegia Cannulae are indicated for the infusion of blood or cardioplegia solution into the coronary venous system, intraoperatively. They are available with either a self- or manually-inflating balloon

| | |
|-----------------------|-----------|
| Order Code | RCM-14315 |
| Size | 15 Fr |
| Available Coated | no |
| Effective Length (cm) | 30.5 |
| Qty Per Box | 10 |
| Latex-Free | yes |
| Connector | Luer lock |