

# Cardioplegia Cannulae, Retrograde, Self-Inflating PVC Balloon, 14 Fr, RCS-13114



The Retrograde Cardioplegia Cannulae are indicated for the infusion of blood or cardioplegia solution into the coronary venous system, intraoperatively. They are available with either a self- or manually-inflating balloon

<b>ORDER CODE</b>	RCS-13114
<b>SIZE</b>	14 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	30.5
<b>QTY PER BOX</b>	5
<b>LATEX-FREE</b>	yes
<b>CONNECTOR</b>	Luer lock