


# Cardioplegia Cannulae, Vessels, VC-11110



The Vessel Cannula is indicated for use in delivery of cardioplegia solution during cardiopulmonary bypass surgery, or to help check for leaks in a harvested vein which will be used for a graft. The Vessel Cannula is a flexible, tapered, polyvinyl chloride tube with a soft distal tip. The distal tip has two sizes of anchor hubs (3mm distal - 3.6mm proximal) designed to accommodate vessels of various diameters. The proximal female luer connection is designed to securely connect with a standard cardioplegia delivery line.

<b>ORDER CODE</b>	VC-11110
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	6
<b>QTY PER BOX</b>	30
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes