

# MICS Cannulae, Aortic Root Long Needle, 12 Ga (9 Fr), AR-17012



Livanova Aortic Root Long Needle is specifically designed for minimally invasive cardiac surgery in lateral thoracotomy, especially mitral valve procedures. It's longer than a conventional aortic root (35 cm) which allows it to be positioned in lateral thoracotomy.

Order Code	AR-17012
Size	12 Ga (9 Fr)
Insertion Needle	14 Ga - Green
Available Coated	no
Effective Length (cm)	31.8
Qty Per Box	10
Latex-Free	yes