

# MICS Cannulae, Aortic Root Long Needle, 14 Ga (7 Fr), AR-17014



Livanova Aortic Root Long Needle is specifically designed for minimally invasive cardiac surgery in lateral thoracotomy, especially mitral valve procedures. It's longer than a conventional aortic root (35 cm) which allows it to be positioned in lateral thoracotomy.

Order Code	AR-17014
Size	14 Ga (7 Fr)
Insertion Needle	16 Ga- White
Available Coated	no
Effective Length (cm)	31.8
Qty Per Box	10
Latex-Free	yes