

# MICS Cannulae, Easy Flow DUO Femoral Arterial, 23 Fr, 103-310



Livanova EasyFlow DUO Femoral Arterial Cannulae feature a dispersion tip for gentle flow in the femoral artery. The obturator allows bloodless insertion in the arterial vessel. Available in the femoral configuration with a dedicated kit for Seldinger insertion. Validated also for venous drainage support.

ORDER CODE	103-310
SIZE	23 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	28,5
QTY PER BOX	5
LATEX-FREE	yes



**Contact us** Complete the form to get in touch with a LivaNova representative.

Contact Us