

# MICS Cannulae, RAP Femoral Venous, 23 Fr Distal/25 Fr Proximal, 200-150



Livanova RAP Femoral Venous Cannulae are designed for improved venous drainage, with a dual stage tip which drains from superior and inferior vena cava. Available in two sizes 22/22 Fr and 23/25 Fr with full-flow drainage capability. A high-quality soft introducer features a smooth transition of the cannula to improve insertion and reduce vascular damage.

|                              |                             |
|------------------------------|-----------------------------|
| <b>ORDER CODE</b>            | 200-150                     |
| <b>SIZE</b>                  | 23 Fr Distal/25 Fr Proximal |
| <b>AVAILABLE COATED</b>      | no                          |
| <b>EFFECTIVE LENGTH (CM)</b> | 65,5                        |
| <b>QTY PER BOX</b>           | 5                           |
| <b>LATEX-FREE</b>            | yes                         |