


MICS Cannulae, Stand-Alone Guidewire, 0.035 (in), IQ35F180J6F



Livanova Vascular Dilator Kit, for percutaneous cannulation technique, includes an 18 Gauge needle, a 0.035 x 180 cm guide-wire, a blade scalpel and 8,12,16,20 and 24 Fr dilators.

ORDER CODE	IQ35F180J6F
SIZE	0.035 (in)
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	180
QTY PER BOX	10
LATEX-FREE	yes



**Contact us** Complete the form to get in touch with a LivaNova representative.

Contact Us