

Suction Products, Suction Wands, Rigid Fluted Tip with Metal Tube and Plastic Handle, 21.7 Fr, SU-12202



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

| | |
|-------------------------|----------|
| ORDER CODE | SU-12202 |
| SIZE | 21.7 Fr |
| AVAILABLE COATED | no |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |