

## Suction Products, Suction Wands, Rigid Fluted Tip with Plastic Handle, 21 Fr, S351-70



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

|                         |                                   |
|-------------------------|-----------------------------------|
| <b>ORDER CODE</b>       | S351-70                           |
| <b>SIZE</b>             | 21 Fr                             |
| <b>AVAILABLE COATED</b> | Depending on country registration |
| <b>QTY PER BOX</b>      | 10                                |
| <b>LATEX-FREE</b>       | yes                               |