

Suction Products, Sumps, ATRASump, 1/4", SU-29602



Livanova Suckers and Sumps are intended for use to remove excess fluid from the surgical field. They are available with flexible or malleable tubing, with weighted metal spring wound tip that terminates in a rounded metal tip. ATRASump offers a unique, a-traumatic, weighted tip with a clear cover available over the metal coil, designed with safety in mind

ORDER CODE	SU-29602
SIZE	1/4 in
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	38
QTY PER BOX	10
LATEX-FREE	yes