

## Thoracic Catheters, Right Angle, 12 Fr, TC-11112



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2" increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

<b>ORDER CODE</b>	TC-11112
<b>SIZE</b>	12 Fr
<b>EFFECTIVE LENGTH (CM)</b>	51,8
<b>LATEX-FREE</b>	yes
<b>CONNECTOR ACCEPTANCE</b>	3/16 in