

Thoracic Catheters, Right Angle, 40 Fr, TC-11140



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2" increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

| | |
|------------------------------|----------|
| ORDER CODE | TC-11140 |
| SIZE | 40 Fr |
| EFFECTIVE LENGTH (CM) | 51,8 |
| LATEX-FREE | yes |
| CONNECTOR ACCEPTANCE | 1/2 in |