

Thoracic Catheters, Straight, 20 Fr, TC-10120



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2''' increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

| | |
|-----------------------|----------|
| ORDER CODE | TC-10120 |
| SIZE | 20 Fr |
| EFFECTIVE LENGTH (CM) | 51,8 |
| LATEX-FREE | yes |
| CONNECTOR ACCEPTANCE | 1/4 in |



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us