

Thoracic Catheters, Straight, 40 Fr, TC-10140



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2''' increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

ORDER CODE	TC-10140
SIZE	40 Fr
EFFECTIVE LENGTH (CM)	51,8
LATEX-FREE	yes
CONNECTOR ACCEPTANCE	1/2 in



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us