

# Venous Return Cannulae, Dual Stage, Low Profile, Wire-reinforced Tubing, 36/50 Fr, LRD-61050



The Dual Stage Venous Return Cannula is indicated for single tube venous drainage from the right atrium and vena cava during cardiopulmonary bypass surgery. The Dual Stage Cannula has a proximal flexible polyvinyl chloride tubing segment for drainage of the right atrium. The catheter body tapers to a flexible polyvinyl chloride tubing distal segment that ends in a straight fluted distal tip for drainage of the vena cava. The proximal and distal segments of the catheter are available in wire or non-wire reinforced tubing.

ORDER CODE	LRD-61050
SIZE	36/50 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	41
QTY PER BOX	10
LATEX-FREE	yes



**Contact us** Complete the form to get in touch with a LivaNova representative.

Contact Us